Cathejell Lidocain C

Sterile lubricant gel

Product information leaflet

en

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cs

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hu

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Návod na použitie

sk

Pharm. Fabrik Montavit Ges.m.b.H.
Salzbergstr. 96, A-6067 Absam, Austria

+43 5223 57926 0

+43 5223 52294

@ pharma@montavit.com

Rev. April 2020

CE 0124

Cathejell Lidocain C

Sterile lubricant gel

Cathejell Lidocain C is a sterile, water-soluble, visually clear gel, which is used as a lubricant for catheters, endoscopes, probes, tubes or other medical instruments.

Cathejell Lidocain C is used

- for instillation into the urethra prior to catheterisation and the introduction of instruments
- use as a lubricant when introducing probes or endoscopes and in tracheal intubation (use in procedures in the respiratory tract, change of suprapubic catheter, in endoscopic procedures in the rectum and colon, as well as in endoscopic applications in the oro- and nasopharyngeal area and in gynaecological applications)

Cathejell Lidocain C serves as a lubricant for use on mucous membranes and, in addition to its principal intended action, has a local anaesthetic and germicidal action; it is intended to relieve pain during such procedures and to prevent infections to a large extent. The effect starts as early as 5 – 10 minutes after application and lasts for 20 – 30 minutes.

During urethral catheterisation, instillation with *Cathejell Lidocain C* gently opens the urethra prior to inserting a catheter or other instruments and forms a lubricant film between the urethra mucosa and the instrument. Thus, the lubricant gel passes through the urethra ahead of the instrument and thereby makes it easier to insert instruments without pain.

During catheterisation, the patient relaxes the pelvic floor muscles, thereby allowing the instrument to pass more easily from the urethra through the external sphincter.

In order to prevent infections within the area of application, *Cathejell Lidocain C* contains chlorhexidine dihydrochloride. This agent disinfects and acts against a large number of bacteria and some fungi and viruses.

Directions for use:

Use in the urethra before inserting a catheter, endoscope or other medical instruments

Important: Slowly apply (instil) into the urethra before inserting any instruments.

The accordion syringes contain 12.5 g or 8.5 g gel, of which approximately 10 g and 6 g, respectively, are inserted into the urethra during instillation.

Method of administration:

- Clean or disinfect, if appropriate, the external opening of the urethra
- Peel off the paper from the transparent cover up to the tapered end of the blister
- Snap off the tip with a short sharp jabbing action preferably into the blister pack
- Remove the tip completely so that it cannot be inserted accidentally into the urethra
- Squeeze out one drop of gel so that the application nozzle can be inserted more easily
- Slowly insert the gel by applying moderate pressure on the syringe.

Use with endoscopes and probes

The accordion syringes contain 12.5 g or 8.5 g gel, of which the designated amount of up to about 10 g and 6 g, respectively, is applied.

- Peel off the paper from the transparent cover up to the tapered end of the blister
- Snap off the tip with a short sharp jabbing action preferably into the blister pack
- Remove the tip completely so that it cannot be inserted by mistake
- Slowly apply the gel by exerting moderate pressure on the syringe.

Apply gel to the probe / endoscope / instrument depending on the intended use and spread evenly.

Use in Airway management

Use in tracheal intubation

The accordion syringes contain 12.5 g or 8.5 g gel, of which the designated amount of up to about 10 g and 6 g, respectively, is applied.

- Peel off the paper from the transparent cover up to the tapered end of the blister
- Snap off the tip with a short sharp jabbing action preferably into the blister pack
- Remove the tip completely so that it cannot be inserted by mistake
- Slowly apply the gel by exerting moderate pressure on the syringe.

Spread the gel (about 5 g) evenly over the lower third of the tube to pass smoothly into the trachea. To reduce the risk of microaspiration, lubricate the cuff prior to intubation. The gel can also be applied on the introducer to facilitate a better sliding into and out of the tube. To prevent drying out, the gel should not be applied to the instrument until immediately before use. Do not let any gel enter the lumen of the tube.

The syringes are intended for single use; discard any remaining gel.

Use with laryngeal mask airway (LMA)

The accordion syringes contain 12.5 g or 8.5 g gel, of which the designated amount of up to about 10 g and 6 g, respectively, is applied.

Method of administration:

- Peel off the paper from the transparent cover up to the tapered end of the blister
- Snap off the tip with a short sharp jabbing action preferably into the blister pack
- Remove the tip completely so that it cannot be inserted by mistake
- Slowly apply the gel by exerting moderate pressure on the syringe.

For correct application spread a thin layer of the gel on the posterior surface of the laryngeal mask. To reduce the risks of blockage of the aperture or inhalation of lubricant, avoid remaining globules of gel on the anterior surface of the cuff or in the bowl of the mask.

The syringes are intended for single use; discard any remaining gel.

Dosage instructions:

The syringes contain 12.5 g or 8.5 g gel, of which about 10 g and 6 g, respectively, can be removed.

The amount of gel to be used depends on the patient's anatomical characteristics and the type of endoscope, catheter, tube, probe or instrument inserted.

For catheterisation, never instil more than 1 syringe at any one time. The contents of one syringe are sufficient to fill the urethra.

In children below 12 years of age, doses should not exceed 2.9 mg lidocaine (=0.14 ml gel) / kg BW. No more than four doses should be administered within 24 hours.

In children below 2 years of age, *Cathejell Lidocain C* must not be used.

In patients at risk, elderly, debilitated and acutely ill patients and in cases of impaired liver function or severe kidney dysfunction, the dosage should be adjusted accordingly. The maximum dose must be calculated in mg lidocaine hydrochloride/kg body weight (2.9 mg lidocaine hydrochloride/kg BW).

Do not use *Cathejell Lidocain C*

- if you are allergic (hypersensitive) to the active substances or any of the other ingredients of *Cathejell Lidocain C*,
- if you are allergic (hypersensitive) to certain other local anaesthetics (of the amide type),
- in children below 2 years of age,
- in bulbocavernosus reflux (an injury to the thin urethral lining, which may cause the lubricant to infiltrate into erectile tissue, where absorption may occur),
- in patients with a severely weak heart, extremely slow heartbeat, heart conduction disturbances (AV block), shock due to heart failure or reduced blood volume.

Take special care with *Cathejell Lidocain C*

- if you have inflammation, injury or an ulcer of the urethra
- if your liver or kidney function is severely impaired
- if you have heart or respiratory tract dysfunction
- if you are elderly, in a weak state of health or acutely ill
- if you are prone to seizures or epilepsy
- if you are suffering from a certain muscle disease (myasthenia gravis)
- if you are on treatment with certain medicines for heart rhythm disorders, known as class III antiarrhythmics (e.g. amiodarone), as the effects on the heart may be enhanced
- in children, as absorption of lidocaine into the blood stream may be increased.
- in patients with acute porphyria or methaemoglobinaemia. In this case, you must seek advice from a doctor before use.

Swallowing difficulties may occur when using *Cathejell Lidocain C* in the mouth and throat region. There is a risk of inhaling the gel (aspiration). Numbness of the tongue or mouth area can lead to a bite injury. Therefore food should be avoided for three hours after the application of *Cathejell Lidocain C* in the mouth and throat region.

If the contents of more than one syringe are inserted into the urethra, if large amounts of gel enter the bladder or if the urethra is ulcerous / inflamed, this may generally cause increased absorption of lidocaine and hence lead to an overdose, particularly in children and elderly patients, with subsequent disorders of the central nervous and cardiovascular system.

Cathejell Lidocain C must not come into contact with the eyes, middle ear or wounds.

In general anaesthesia, a lubricant without lidocaine should be preferred.

The use of anionic tensides (certain soaps) can considerably reduce the disinfectant action of chlorhexidine.

Interaction between *Cathejell Lidocain C* and other medicines

Due to possible increased effects on the heart, *Cathejell Lidocain C* should not be used at the same time as:

Medicines containing lidocaine or other amide-type local anaesthetics, antiarrhythmics (used to treat heart rhythm disorders) or beta-blockers or calcium channel blockers (used for high blood pressure).

Taking cimetidine (inhibits stomach acid production) can increase the risk of side effects.

Pregnancy and breast-feeding:

During pregnancy and breast-feeding, *Cathejell Lidocain C* should be used only after careful consideration of the benefits and risks by a doctor and after the individual dose has been established.

Repeated use during pregnancy and breast-feeding is not recommended.

Please note that there should be a 12-hour interval between using *Cathejell Lidocain C* and subsequent breast-feeding.

Women planning a pregnancy or who think they are pregnant should tell their doctor.

Driving and using machines

Effects on the ability to drive and use machines are unlikely, but cannot be completely ruled out in cases of increased individual sensitivity.

Side effects:

Side effects rarely occur after *Cathejell Lidocain C*, provided that the product is used according to the dosage recommendations/recommendations for use and the necessary precautions are taken.

In rare cases, allergic reactions such as redness, burning or itching may occur. However, more severe, general hypersensitivity reactions and even shock are also possible.

In view of the low exposure levels and short duration of use with *Cathejell Lidocain C*, other significant side effects/interactions are unlikely.

Systemic side effects may be caused by high plasma levels, rapid absorption or an overdose, or by hypersensitivity, idiosyncrasy or reduced tolerance, whereupon the following symptoms may occur:

CNS: nervousness, dizziness, blurred vision, tremor. These signs may not necessarily occur; in some patients, intoxication manifests as drowsiness, unconsciousness or respiratory arrest.

Cardiovascular: hypotension, bradycardia, asystole.

As a local reaction, hoarseness can occur during use as a lubricant for the endotracheal tube.

Special warnings for safe use:

In anaesthesia, a lubricant without lidocaine should be preferred.

Only use packs with the sterile paper intact.

After instillation, a 5 to 10-minute interval must be respected before pain-free insertion of the instrument is possible.

If reused or resterilised, sterility, functionality and product quality are no longer assured.

Composition:

1 g contains 20 mg lidocaine hydrochloride and 0.5 mg chlorhexidine dihydrochloride.

The other ingredients are: hydroxyethylcellulose, glycerol, sodium lactate, water for injections.

Pack sizes:

1 x 12.5 g, 5 x 12.5 g, 25 x 12.5 g

1 x 8.5 g, 5 x 8.5 g, 25 x 8.5 g

Shelf life:

As printed on the pack.

Storage instructions:

Do not store above 25°C.

Keep the blister in the carton in order to protect from light.

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